

## US EPA approved Microbiological Methods for Testing Drinking, Ground, Recreational and Waste Waters

Gil Dichter  
World Wide Technical Support Manager, Water  
[www.idexx.com/water](http://www.idexx.com/water)



---

---

---

---

---

---

---

---

## OBJECTIVES

- Introduction to water microbiology
- Define total coliforms, fecal, E.coli & enterococci
- Overview of Regulations
- USEPA approved methods
- Understanding MPN theory and relationship to MF
- Q&A



---

---

---

---

---

---

---

---

## Bacteria Hierarchy

- Family
- Genus (type of bacteria)
- Species
- Enterobacteriaceae
- Escherichia, Klebsiella, Citrobacter
- E.coli, E. hermanii, K. pneumoniae, C. freundii



---

---

---

---

---

---


---

---

### Coliform Genera

- Escherichia- human and animal feces
- Enterobacter- environment, feces
- Klebsiella- environment, feces
- Citrobacter- environment
- Serratia- environment

---




---

---

---

---

---

---


---

---

### Thermotolerant Coliforms [aka Fecal Coliforms]

- It is a subset of total coliforms
- Defined as coliform bacteria that can grow at 44.5°C
- Consists of the following coliforms:
  - *E. coli*
  - *K. pneumoniae*
  - *Enterobacter spp.*

---




---

---

---

---


---

---

---


---

### Coliform Bacteria Group



The diagram illustrates the relationship between different groups of coliform bacteria. A large green circle represents 'Total Coliforms (35°C)'. Inside it, a smaller blue circle represents 'Fecal Coliforms-thermotolerant (44.5°C)'. Within the blue circle, there are two smaller white circles: 'Klebsiella' on the left and 'E.coli' on the right. Below the blue circle, the words 'Citrobacter' and 'Enterobacter' are listed, indicating they are part of the total coliforms but not necessarily part of the fecal coliforms group.

---




---

---

---

---

---

---

---

---

## Definitions of Coliforms & E.coli

- Standard Methods for the Examination of Water and Wastewater, 21<sup>st</sup> Edition
  - **9221 MTF:** Gram negative bacteria that ferment lactose resulting in gas and/or acid formation (turbidity) within 48 h at 35C.
  - **9222 MF:** Gram negative bacteria that develop red colonies with a metallic sheen within 24 h at 35C on m-Endo medium. Some members of the coliform group produce red colonies without a metallic sheen.
  - **9223 Enzymatic:** The total coliform group is defined as all bacteria possessing the enzyme β-D galactosidase cleaving the chromogenic substrate resulting in the release of the chromogen.




---

---

---

---

---

---

---

---

## Definition of coliforms and E.coli

Enzymatic con't

*E.Coli* is defined as giving a total coliform response and possessing the enzyme β-D glucuronidase, cleaving the fluorigenic substrate resulting in the release of the fluorgen.

- Fecal Coliform (thermotolerant bacteria)
  - Bacteria that grow at 44.5°C
    - Produce gas & turbidity in EC medium
    - Produce blue colonies in m-FC media
    - Produce yellow color wells with Colilert-18




---

---

---

---

---

---

---

---

## Why Test for E.coli

Animal	# Tested	<i>E.coli</i>	<i>Klebsiella</i> spp	Enterobacter/ Citrobacter
Human	26	96.8%	1.5%	1.7%
Cow	15	99.9	-	0.1
Horse	3	100	-	-
Sheep	20	97	-	3
Pig	15	83.5	6.8	9.7
Average		94.5%		

Source: *E.coli*: Fecal Coliform A.P. Dufour. Special Technical Publication 65, ASKIM, Pp48-58, 1977




---

---

---

---

---

---

---

---

### Commonly Used Indicator Bacteria for Water Testing

- Most commonly used:
  - Total coliform
  - Fecal coliform
  - *E. coli*
  - *Enterococci*
- Total coliform used for >100 years
- Fecal coliform used for >80 years
- *E. coli* >20 years

10



---

---

---

---

---

---

---

---

### Requirements for an Indicator Organism

- Present when pathogens are present
- Absent in uncontaminated waters
- Present in higher numbers than pathogens in contaminated water
- Better survival in water than pathogens
- Easy and Safe to analyze
- Rapid results
- Inexpensive
- Accurate

11



---

---

---

---

---

---

---

---

### Overview of Regulations

- Total Coliform Rule-
- LT2
- Ground Water
- Recreational Water
- Waste Water

12



---

---

---

---

---

---


---

---

## Total Coliform Rule

Promulgated 6/29/89

- To control pathogens
- Monitor for total coliform bacteria in the distribution system (CFR 141)
  - Number of samples dependent on population served
    - i.e.; 8501-12900 population requires 10 samples/month minimum.
- Used as an indicator of treatment reliability & distribution integrity
- If >5% TC/month, violates MLC for TC; if 2 consecutive TC with one E.coli, = acute violation and requires immediate notification



---

---

---

---


---

---

---


---

## Where is Sampling Done?




**Large Systems**

- May have dedicated sampling stations.
- In public buildings such as fire and police stations.
- Some are sampled daily M-F while others once a week, or two or three times a week.



**Small Systems**

- Depending on size as little as one site exists for sampling.
- Usually no dedicated sampling stations
- Maybe a kitchen or laundry sink.



---

---

---

---

---


---

---

---

## LT2 ESWTR

- Proposed- 7/03
- Promulgated- 1/5/06
- Large water systems must test for Crypto, E.coli & turbidity
- Small water systems test for E.coli
- Small systems may be required to test for Crypto.
- Proposes methods for Crypto & E.coli



---

---

---

---

---

---

---

---

## Key Features of the LT2 Rule

- Applies to all public water systems using surface water sources (including GWUDI)
- Supplements existing regulations to address *Cryptosporidium* in systems with higher risk
  - Filtered systems with high source water occurrence
  - All unfiltered systems

15



---

---

---

---

---

---

---

---

## Analytical Methods for Biological Pollutants in the Ambient Water

### Rule

- Procedure for enumerating E.coli, enterococci and Crypto & Giardia in ambient water
- Based on 1986 guidelines for E.coli & enterococci
- Methods include- m-TEC, Colilert and Colilert-18 for E.coli, m-EI & Enterolert for enterococci.
- Final rule effective - 8/20/03

16



---

---

---

---

---

---

---

---

## Guidelines for Analysis of Pollutants; Analytical Methods for Biological Waste Water

- Federal Register Notice- August 16<sup>th</sup>, 2005
- Promulgated on March 26<sup>th</sup>, 2007
- Testing for E.coli and for Enterococci
  - Approved methods include m-TEC, m-ColiBlue, Colilert, Colilert-18 for E.coli, Enterolert & m-EI methods for enterococci
  - E.coli for fresh surface waters and enterococci for marine waters
- June 2010**
- US EPA recommends Colilert-18 for testing fecal coliforms in WW AT 45.5°C

17



---

---

---

---

---

---

---

---

## Multiple Tube Fermentation (MTF) Most Probable Number (MPN)

15




---

---

---

---

---

---

---

---

---

---

---

---

### Most Probable Number Assay

**Five-Tube MPN Table**

NUMBER OF TUBES GIVING POSITIVE REACTION	MPN INDEX PER 100 ML		NUMBER OF TUBES GIVING POSITIVE REACTION	MPN INDEX PER 100 ML		
	0	1			0	1
0	0	0	4	2	1	26
0	0	1	4	3	0	22
0	1	0	4	3	1	33
0	1	1	4	4	0	34
1	0	0	3	5	0	23
1	0	1	4	5	0	30
1	1	0	4	5	1	40
1	1	1	4	5	1	30
2	0	0	3	5	1	40
2	0	1	3	5	2	30
2	1	0	3	5	2	40
2	1	1	3	5	2	40
2	2	0	2	5	3	140
2	2	1	3	5	3	140
3	0	0	2	5	3	170
3	1	0	2	5	4	130
3	1	1	2	5	4	170
3	2	0	2	5	4	220
3	2	1	2	5	4	260
4	0	0	1	5	4	290
4	0	1	1	5	5	340
4	1	0	1	5	5	340
4	1	1	1	5	5	340
4	2	0	1	5	5	1,400
4	2	1	1	5	5	1,400

NOTE: 0 = a number of broth tubes showing gas after inoculation with 10 ml sample.  
1 = a number of broth tubes showing gas after inoculation with 1 ml sample.  
2 = a number of broth tubes showing gas after inoculation with 0.1 ml sample.  
SOURCE: Standard Methods for the Examination of Water and Wastewater, 19th ed. American Public Health Association, 1995.

20




---

---

---

---

---

---

---

---

---

---

---

---

## Multiple Tube Fermentation Method (Total Coliforms-LTB & BGLB)



21




---

---

---

---

---

---

---

---


---

---

---

---

**Fecal Coliform/E. coli Multiple Tube Method -  
EC medium + MUG**



21

**IDEXX**

---

---

---

---

---

---

---

---

**Membrane Filtration  
Methods**

22

**IDEXX**

---

---

---

---

---

---

---

---

**Membrane Filter Apparatus**



23

**IDEXX**

---

---

---

---

---

---

---

---



### Membrane Filtration- applicable to all methods

- May require media prep
- Extensive QC
- Minimum of 20 steps
- 20-80 colonies (20-60 for m-FC)
- Risk of confluent growth
- Risk of clogged filters
- Risk of overlapping colonies
- Risk of air bubbles under membrane
- Difficult to read

25



---

---

---

---

---


---

---

---

### m-Endo: Method

- Standard Methods 9222B&C
- Incubate for 24 hours at 35-37°C
- Read at 24 hours for presumptive results
  - Dark red or pink with metallic/green sheen
- HPC Interference at 120 colonies or more
  - Transfer colonies to BGB to confirm coliforms (up to 48 hours)
    - Positives produce gas and turbidity
- Transfer positives to EC (+Mug) for 24 hours at 44.5°C
  - Positives produce, turbidity/gas, and/or fluorescence
- Results with confirmation- could take up to 3 days



26



---

---

---

---


---

---

---

---

### m-Endo Plate



27



---

---

---

---

---


---

---

---

**Hach m-ColiBlue**

- Reagent storage is at 2-8°C
- Incubate for 24 hours at 35±0.5°C
- Read at 24 hours
  - Red or blue colony = TC
  - Blue colony = *E.coli*
  - Difficult to interpret colonies (very small)
- If count >200 for TC/EC, exclude from calculation



**IDEXX**

---

---

---

---

---


---

---

---

**Modified m-TEC for E. coli**

- Contains 5-bromo-6-chloro-3-indolyl-β-D-glucuronide
- Selective inhibitory chemicals that can effect the growth of sub-lethal injured bacteria.
- 2 hours of incubation at 35°C followed by 22 hours at 44.5°C
- Positive reaction is red to magenta color colonies
- Extensive QC requirements



**IDEXX**

---

---

---

---

---


---

---

---

**m-FC**

- Enriched lactose medium
- Contains aniline blue
- Larger colonies than m-Endo; range is 20-60 colonies
- Incubate in water bath at 44.5 ±0.2°C for 24 hours in plastic bag
- Blue colonies are positive
- Atypical- grey to green colonies
- Confirm as per section 9020 Standard Methods



**IDEXX**

---

---

---

---

---

---

---

---

## Defined Substrate Technology- Colilert & Colilert-18

25

IDEXX

---

---

---

---

---


---

---

---

### ONPG Positive Reaction

Colilert & Colilert-18



26

IDEXX

---

---

---

---

---


---

---

---

### MUG Positive Reaction

Colilert & Colilert-18



27

IDEXX

---

---

---

---

---

---

---

---

### MUG Positive Reaction

Colilert & Colilert-18

**E. coli**  
 $\beta$ -glucuronidase  
 $\beta$ -D-glucuronide  
 MUG  
 4-methyl-umbelliferyl

**E. coli**  
 $\beta$ -glucuronidase  
 $\beta$ -d-glucuronide  
 4-methyl-umbelliferone

24

**IDEXX**

---

---

---

---

---

---

---

---

---

---

### Quanti-Tray Demonstration

Add Colilert to sample and shake to dissolve

Pour mixture into a Quanti-Tray

25

**IDEXX**

---

---

---

---

---

---

---

---

---

---

### Quanti-Tray Demonstration

Seal and incubate at 35°C for 18 or 24 hours for total coliforms and 44.5°C for fecal coliforms for 18 hours

Count positive wells and refer to MPN table

26

**IDEXX**

---

---

---

---

---

---

---

---

---

---

**E.coli- Blue Fluorescence- Quanti-Tray**  
under a 365nm UV Light at 35°C for 18 or 24 hours for  
E.coli



27



---

---

---

---

---

---

---

---

---

---

### Studies

➤ AWWARF Study : **Significance of Methods and Sample Volumes for EC & TC Measurements- Fricker et al**

– Methods and volume of samples evaluated for drinking water testing cited in Sept. 2010 RTCR meeting.

- Glucosidase methods will detect more coliforms (higher concentration and a broader range) than lactose based methods
- Study found Colilert-18 and Colilert had a confirmation rate of 97.7 and 95.1%. False negative rate of 1.8 and 2.8%
- Confirmation rates for m-Endo was 84.6% & false negative rate of 13.2%
- Verification of E.coli for Colilert-18 & Colilert are 100 & 99.1%, m-FC is 79.7%

28



---

---

---

---

---

---

---

---

---

---

### Studies con't

➤ Choice of a specific method is often a personal preference within a lab. Factors such as sensitivity, specificity, accuracy, simplicity & time to results are most important factors in choosing an appropriate method

29



---

---

---

---

---

---

---

---

---

---


**Studies con't**

➤ Comparison of m-Endo, MacConkey & Teepol Agar Methods for MF Counting Total Coliform bacteria in Water: Grabow & Preez AEM, Sept'79 Vol 38 #3 p 351-355

- 341 coliform like colonies picked from m-Endo membranes for identification
  - 186 /55% true coliforms
  - 138/40% Aeromonas hydrophilia
  - 17/5% other species

➤ A comparison of ten USEPA approved total coliform/E. coli tests: Jeremy Olstadt, WQTC 2005

- 2 different strains of Aeromonas spike into different waters at 3 different sites- Colilert & Colilert 18 negative at 10<sup>5</sup>/mL



---

---

---

---

---

---

---

---

---

---


**Studies con't**

➤ Evaluation of Two Colilert Formulations and Quanti-Tray for the Assessment of the Bacteriological Quality of Water: Fricker et al. Water Research 1997 Vol 31 pp 2495-2499

- Identified 296 wells for total coliforms for Colilert & Quanti-Tray
  - All 296 wells confirmed as coliforms compared to 257 colonies/363 (71%) for the UK MF method.

➤ Evaluation of the Autoanalysis Colilert for Determination & Enumeration of Coliforms: Covert et al AEM Oct 1989 Vol 55 #10

- No false positives from Pseudomonas, 1 + Aeromonas\* at 7.5X10<sup>4</sup>/ml
- Flavobacterium + at 10<sup>7</sup>-10<sup>9</sup>/mL
- \*All above label claim



---

---

---

---

---

---

---

---

---

---

**Enterococci  
MF and Defined Substrate  
Technology®**



---

---

---

---

---

---

---


---

---

---

### m-EI Method for Enterococci

- 24 hour MF Test
- All colonies regardless of color with a blue halo
- 20-60/plate
- Colonies <0.5mm difficult to read
- Media is expensive (Indoxyl-β-D-glucoside)
- Magnification required to read plates



---

---

---

---


---

---

---

---

### Method 1600- m-EI



---

---

---

---

---

---

---

---

### Enterolert™ Demonstration



Add reagent



Seal in a Quanti-Tray  
Incubate at 41.5C for 24 hrs



---

---

---

---


---

---


---

---

### Enterolert™ Demonstration



Count fluorescent wells  
and refer to MPN table



---

---

---

---

---

---

---

---


---

---

### Estimation of Bacterial Densities by MPN

- Dates back to 1915 when the concept was introduced by M.H. McCrady (J of Infectious Disease Vol 17, 1915\*)
- Prior to this novel concept; no means of direct counting
  - Only had presence-absence of fermentation tubes
- The method is a means for estimating without any direct count, the density of organisms in a liquid.
- Multiple samples of the liquid are taken and incubated in suitable media
  - Record presence or absence of growth in each sample tube
  - Ingenious application of probability theory
  - Estimate the number of organisms from the number of negative tubes

\*The Numerical Interpretation of Fermentation Tube Results



---

---

---

---

---

---

---

---

---

---

### MPN Theory – a short math lesson


Basic assumptions:

- The organisms are distributed randomly throughout
- the liquid
  - Sample is well shaken (important and often neglected)
- A sample will exhibit growth (in the culture media) whenever one or more of the target organisms is present
  - Media should be selective and sensitive
- Requires that at least 1 tube shows no growth (sterile)

If  $n$  samples, each of volume  $v$  mL are taken from a liquid, and  $s$  of these are sterile, then an estimate of the organism density  $d$  per mL in the original sample is:

$d = -2.303/v \log(s/n)$  [Poisson distribution]

$d$  can be shown statistically to have the highest probability of estimating the actual organism density – the "Most Probable Number" - MPN



---

---

---

---

---

---

---

---

---


---





### Overcoming the tube numbers problem

- To increase counting range need more tubes
- For 100mL sample could use 50 tubes x 2mL
- Would give counting range of ~200MPN/100mL
  - Big headache for lab – lots of tubes – lots of pipettes -lots of washing! Not really practical.
  - Microtitre plate (96 wells) – MPN range 438/100mL




---

---

---

---

---

---

---

---

---

---

### Quanti-Tray®

- Use multiple well concept but automate the sealing



**51 well Quanti-Tray**  
Total sample volume 100mL

$$d = -2.303/v \log(s/n)$$

$$n = 51, v = 1.96 \text{ mL}, s = 1$$

$$d = -2.303/1.96 \log(1/51) \\ = 2.006 \text{ MPN/mL} \\ = 200.6 \text{ MPN/100 mL}$$




---

---

---

---

---

---

---

---

---

---

### Multiple Dilutions

- Not practical to extend the counting range by increasing the number of wells beyond a certain point.
- Sample serially diluted and each dilution is inoculated into a similar number of tubes. Commonly used configurations:
  - 3 tubes x 3 decimal dilutions – food industry
  - 5 tubes x 3 decimal dilutions – food, water wastewater (Max MPN = 1600/100mL)
- With multiple dilutions the statistical calculation of the MPN value becomes more complex. **Thomas** provided an **approximation** that could be used for any combination of tubes

$$\text{MPN/100mL} = \frac{\# \text{ of positive tubes} \times 100}{\sqrt{(\# \text{ of mL in negative tubes}) \times (\# \text{ of mL in all tubes})}}$$




---

---

---

---

---

---

---

---

---

---



### 95% Confidence Limit Comparison

Quanti-Tray MPN/100mL	Quanti-Tray 2000 MPN/100mL	MTF MPN/100mL	MF CFU/100mL
Value CL	Value CL	Value CL	Value CL
5.3 2.3-12.3	5.2 1.8-10.8	4.5 0.79-15	5 1.6-7.2
17.8 10.8-29.4	17.9 10.7-28.2	17 6.0-40	17 9.9-27.2
20.7 13-33.3	20.6 12.7-31.8	20 6.8-40	20 12.2-30.8
50.4 35.4-72.5	50.4 35.0-69.1	49 15-150	
78.2 56.4-111.2	78.0 55.6-103.8	79 22-220	




---

---

---

---

---

---

---

---

### IDEXX Support

1-800-321-0207

- #1 Customer Support
- #2 Technical Service
- #3 Select extension

- [www.idexx.com/water](http://www.idexx.com/water)




---

---

---

---

---

---

---

---

Thank You

Questions




---

---

---

---

---

---

---

---